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



1.

2.

3.

4.

The screenshot shows the HAYEAR software interface. The main window displays a video frame with a ruler at the top and a scale on the left. A calibration dialog box is open in the center, with the following fields:

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Below the fields is a diagram illustrating the calibration process, showing a ruler with a segment labeled "像素" (Pixel) and a longer segment labeled "实际长度" (Actual Length).

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5.

The screenshot displays the HAYEAR software interface. The main window shows a video frame with a ruler at the top and a vertical scale on the left. A calibration dialog box is open in the center, titled "定标" (Calibration). The dialog box contains the following fields and values:

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Below the fields is a diagram showing a ruler with "像素" (Pixel) and "实际长度" (Actual Length) labels. The diagram illustrates the relationship between the pixel count and the actual length.

Annotations in the image include:

- A pink arrow pointing to the vertical scale with the text "示范: 1毫米工作距离" (Demonstration: 1mm working distance).
- A red circle around the value "218.209 像素, 1.000 mm" in the dialog box.
- A pink arrow pointing to the "实际长度" field with the text "使用鼠标拖动, 调节距离" (Use mouse to drag, adjust distance).

The software interface also shows a menu bar with options like "文件(F)", "编辑(E)", "查看(V)", "浏览(B)", "设置(S)", "捕获(C)", "图像(I)", "处理(P)", "图层(L)", "测量(M)", "选项(O)", "窗口(W)", and "帮助(H)". The toolbar includes various tools for capture, zoom, and measurement. The left sidebar contains a "相机" (Camera) section with "相机列表" (Camera List) and "捕获与分辨率" (Capture and Resolution) settings. The bottom left corner has a "测量表格" (Measurement Table) button.

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